


Binding Kinetics of a GPCR

Thomas Raschle[#], Philipp Heine^{*}, Andreas Plückthun^{*}, Matyas Vegh[#] and Kaspar Cottier[#].

[#]Creoptix AG; ^{*}Dep. of Biochemistry, University of Zurich, Winterthurerstr. 190, 8057 Zurich, Switzerland



G-protein coupled receptors (GPCRs) comprise a large class of integral membrane proteins, crucial for the transduction of extracellular stimuli across the plasma membrane to elicit molecular and cellular responses. They regulate a variety of physiological processes in eukaryotes and represent the largest group of therapeutic targets, with more than 30% of available pharmaceuticals targeting GPCRs.

Membrane proteins are notoriously difficult to study due to the requirement for a membrane-mimicking environment and their instability once extracted from a cellular membrane. Here we show the capability of the WAVE to measure the interaction of a peptide ligand agonist (NTA11) with a thermostabilized variant of the neurotensin receptor 1 (NTSR1)¹ at highest resolution. This GPCR mediates the multiple functions of neurotensin, such as hypotension, hyperglycemia, hypothermia, antinociception, and regulation of intestinal motility and secretion².

Material and Methods

The receptor was in vivo biotinylated via an avi-tag and captured on a streptavidin pre-coated sensor (WAVEchip DXH-S). For the measured peptide agonist, a mutated and truncated form of neurotensin (Mw 725 Da) comprising residues 8-13 with a Y11A substitution was used. Dose-response curves were recorded for 7 different analyte concentrations of a 3-fold dilution series with 300 nM being the highest concentration. The flow rate was set to 30 $\mu\text{L}/\text{min}$. The running buffer was 50 mM Tris pH 7.5, 150 mM NaCl, and 0.1% (w/v) of the detergent lauryl maltose-neopentyl glycol (L-MNG)³. All measurements were carried out at 25°C.




Figure 1: Sensorgrams of the interaction between a modified neurotensin peptide (Mw 725 Da) as analyte and a thermostabilized neurotensin receptor (NTSR1, 70 kDa protein/detergent complex) as ligand. The data were fitted with a model for a 1:1 interaction including a term for mass transport limitation in the WAVEcontrol software.

Kinetic determination				Equilibrium determination	
R_{\max} (pg/mm ²)	k_{on} (M ⁻¹ s ⁻¹)	k_{off} (s ⁻¹)	K_d (nM)	R_{\max} (pg/mm ²)	K_d (nM)
17.982	1.681E6	3.497E-2	20.8	22.083	42.6

Table 1: Kinetic rates and dissociation constant for the NTA11 agonist obtained with a 1:1 interaction model including a term for mass transport limitation.

Results

Purified and biotinylated neurotensin receptor was captured on a streptavidin pre-coated sensor (WAVEchip DXH-S) at a density of 1350 pg/mm². Dose-response curves for binding of peptide agonist were recorded (Figure 1), yielding binding data that could be well fitted with a model for a 1:1 interaction including a term for mass transport. The obtained kinetic rates and equilibrium constant are summarized in Table 1.

References

- (1) Egloff, P., Hillenbrand, M., Klenk, C., Batyuk, A., Heine, P., Balada, S., Schlinkmann, K. M., Scott, D. J., Schütz, M., and Plückthun, A. (2014) Structure of signaling-competent neurotensin receptor 1 obtained by directed evolution in Escherichia coli. Proc. Natl. Acad. Sci. U. S. A. 111, E655–62.
- (2) Krumm, B. E., and Grisshammer, R. (2015) Peptide ligand recognition by G protein-coupled receptors. Front. Pharmacol.
- (3) Chae, P. S., Rasmussen, S. G. F., Rana, R. R., Gotfryd, K., Chandra, R., Goren, M. a, Kruse, A. C., Nurva, S., Loland, C. J., Pierre, Y., Drew, D., Popot, J., Picot, D., Fox, B. G., Guan, L., Gether, U., Byrne, B., Kobilka, B., and Gellman, S. H. (2010) Maltose-neopentyl glycol (MNG) amphiphiles for solubilization, stabilization and crystallization of membrane proteins. Nat. Methods 7, 1003–1008.