

SPOTTING THE WEAKEST BINDERS

Fabio Andres¹, Jealemy Galindo¹, Natalia Matassova², Marten Beeg³, Fabio M. Spiga¹

¹ Creoptix AG, Zugerstrasse 76, 8820 Wädenswil, Switzerland, info@creoptix.com

² Vernalis Ltd, Granta Park, Abingdon, CB21 6GB, Cambridge, U.K.

³ Istituto di Ricerche Farmacologiche Mario Negri IRCCS, 20156 Milano, Italy, marten.beeg@marionegri.it

A NEW LEVEL OF PERFORMANCE AND FLEXIBILITY IN DRUG DISCOVERY

Employing our Grating-Coupled Interferometry (GCI)¹ technology to deliver superior sensitivity over traditional Surface Plasmon Resonance (SPR) technologies, researchers can reliably determine off-rates of up to 10s^{-1} , starting with just a crude reaction mixture. Combined with microfluidics that sustain a wide variety of solvents - including acetonitrile and high concentrations of DMSO -, the occurrence of false positives is minimized.

Screen, rank and characterize weak binders with off-rates up to 10s^{-1}

Study binding kinetics even at large analyte:ligand MW ratios (up to 1:1000)

Experiment with crude mixtures, detergents and other additives without clogging

GRATING-COUPLED INTERFEROMETRY (GCI)

GCI is a surface-based, label-free biosensing technique. When target molecules (e.g. proteins) are attached to the sensor surface, binding of analytes leads to an increase in mass and hence to a change in the refractive index within the evanescent field near the surface. In GCI, refractive index changes on a sensor surface are measured as time-dependent phase-shift signals. The long light-to-sample interaction length of the waveguide provides intrinsically high signal-to-noise levels for improved sensitivity.

With high sensitivity, the ability to resolve extremely rapid dissociating kinetics and innate compatibility with high molecular weight ratios, the Creoptix® WAVEsystem's GCI technology improves fragment-based screening and kinetic analysis of small molecules to accelerate drug development. Paired with no-clog WAVEchips®, a wide range of molecules can be immobilized using various chemistries.²




More information about GCI


HITS WON'T BE MISSED

The Creoptix® WAVEsystem provides a versatile platform for fragment screening and analysis. Thanks to its ultra-fast transition times, excellent resolution for compounds with k_{off} 's well above 1s^{-1} can be achieved.

Sensorgrams of selected compounds of a 83 crude non-purified compounds^{3,4} library-subset, screened against His-tagged Pyruvate Dehydrogenase Kinase (PDK-1) and captured on a PCH WAVEchip. Samples were injected at ~20 μM . His-tagged HSP90, captured on a different channel, served as negative control.




Selected HSP90 or PDK-1 hits (highlighted) are in good agreement with the reported literature^{3,4}. Aberrant samples can be simply and automatically excluded based on the apparent k_{off} and the χ^2 (grey triangles).



WEAK BINDERS, STRONG DATA


Weak binders such as those found in fragment-based screening are typically ranked by affinity rather than kinetics due to their very fast off-rates, which can not be resolved by traditional SPR instrumentation. Here we show that the Creoptix® WAVEsystem provides an outstanding resolution whereby very fast kinetics can be reliably determined at off-rates up to 10s^{-1} .




Self-assembled amyloid fibrils were immobilized via amine coupling on a 4PCZ WAVEchip® (zwitterionic surface). The small molecule thioflavin (THT, 319 Da) was injected in four (4) concentrations (50 μM - 6.25 μM) for 30s at 400 $\mu\text{l}/\text{min}$. Raw data were double referenced and globally fit with a 1:1 binding model showing accurate determination of an off-rate around 10s^{-1} .

SMALL MOLECULES CAN'T HIDE ANYMORE

Sensitivity is key and often limiting for accurate and reliable analysis of molecular interactions. The high-sensitivity of the Creoptix® WAVEsystem allows researchers to confidently analyze binding interactions at very low signal levels and high analyte-to-ligand molecular weight (MW) ratios.



Sensorgrams of acetazolamide (222 Da) binding to Carbonic Anhydrase II (29 kDa) immobilised at low density onto a PCH WAVEchip.



Sensorgrams of a small drug molecule (295 Da) binding to a target protein (110 kDa) immobilised at low density onto a PCH WAVEchip. Note the analyte-to-ligand MW ratio is > 300.

GCI IS FEATURED IN

NATURE Mechanistic basis for the activation of plant membrane receptor kinases by SERK-family coreceptors

PLANT PHYSIOLOGY The SERK3 elongated allele defines a role for BIR ectodomains in brassinosteroid signalling

eLIFE Dynamics of human protein kinase Aurora A linked to drug selectivity

REFERENCES

19. Kozma et al., "Grating coupled optical waveguide interferometer for label-free biosensing", Sensors and Actuators B: Chemical, 155:444-450 (2011)
2. Höhmann et al., "The SERK3 elongated allele defines a role for BIR ectodomains in brassinosteroid signalling", Nature Plants, 4:345-351 (2018)
3. Brough et al., "Application of Off-Rate Screening in the Identification of Novel Pan-isoform Inhibitors of Pyruvate Dehydrogenase Kinase", J. Med. Chem., 60 (5), pp 2271-2285 (2017)
4. Brough et al., "4,5-Diaryloxazole Hsp90 Chaperone Inhibitors: Potential Therapeutic Agents for the Treatment of Cancer", J. Med. Chem., 51 (2), pp 190-218 (2008)

Scan this code to request a copy of this poster

