

# SPOTTING THE WEAKEST BINDERS



Maria Walton<sup>1</sup>, Jealemy Galindo<sup>1</sup>, Natalia Matassova<sup>2</sup>, Marten Beeg<sup>3</sup>, Fabio Andres<sup>1</sup>, Fabio M. Spiga<sup>1</sup>

<sup>1</sup> Creoptix AG, Zugerstrasse 76, 8820 Wädenswil, Switzerland, info@creoptix.com

<sup>2</sup> Vernalis Ltd, Granta Park, Abingdon, CB21 6GB, Cambridge, U.K.

<sup>3</sup> Istituto di Ricerche Farmacologiche Mario Negri IRCCS, 20156 Milano, Italy

## A NEW LEVEL OF PERFORMANCE AND FLEXIBILITY IN DRUG DISCOVERY

By bringing together modern label-free technology, application development know-how and sophisticated software, Creoptix offers a unique optical biosensor tool for binding kinetics. Engineered around our proprietary Grating-Coupled Interferometry (GCI)<sup>1</sup> technology, the Creoptix WAVEsystem delivers high-quality kinetic data across a broader range of samples than traditional SPR equipment.

Screen, rank and characterize weak binders with off-rates up to  $10 \text{ s}^{-1}$

Study binding kinetics even at large analyte:ligand MW ratios (up to 1:1000)

Experiment with crude mixtures, detergents and other additives without clogging

## GRATING-COUPLED INTERFEROMETRY (GCI)



GCI is a surface-based, label-free biosensing technique. When target molecules (e.g. proteins) are attached to the sensor surface, binding of analytes leads to an increase in mass and hence to a change in the refractive index within the evanescent field near the surface. In GCI, refractive index changes on a sensor surface are measured as time-dependent phase-shift signals. The long light-to-sample interaction length of the waveguide provides intrinsically high signal-to-noise levels for improved sensitivity.


The Creoptix<sup>®</sup> WAVEsystem combines GCI with innovative no-clog microfluidics, allowing the study of interactions even between very small ligands and large complexes. The system uses a robust microfluidic sensor, the WAVEchip<sup>®</sup>, where (membrane) proteins,<sup>2</sup> antibodies, VLPs, peptides or other molecules can be immobilized using various chemistries.




### HITS WON'T BE MISSED

The Creoptix<sup>®</sup> WAVEsystem provides a versatile platform for fragment screening and analysis. Thanks to its ultra-fast transition times, excellent resolution for compounds with  $k_{\text{off}}$ 's well above  $1 \text{ s}^{-1}$  can be achieved.

Sensorgrams of selected compounds of a 83 crude non-purified compounds<sup>3,4</sup> library-subset, screened against His-tagged Pyruvate Dehydrogenase Kinase (PDK-1) and captured on a PCH WAVEchip. Samples were injected at  $\sim 20 \mu\text{M}$ . His-tagged HSP90, captured on a different channel, served as negative control.




Selected HSP90 or PDK-1 hits (highlighted) are in good agreement with the reported literature<sup>3,4</sup>. Aberrant samples can be simply and automatically excluded based on the apparent  $k_{\text{off}}$  and the  $\chi^2$  (grey triangles).



### WEAK BINDERS, STRONG DATA

Weak binders such as those found in fragment-based screening are typically ranked by affinity rather than kinetics due to their very fast off-rates, which can not be resolved by traditional SPR instrumentation. Here we show that the Creoptix<sup>®</sup> WAVEsystem provides an outstanding resolution whereby very fast kinetics can be reliably determined at off-rates up to  $10 \text{ s}^{-1}$ .




Self-assembled amyloid fibrils were immobilized via amine coupling on a 4PCZ WAVEchip<sup>®</sup> (zwitterionic surface). The small molecule thioflavin (THT, 319 Da) was injected in four (4) concentrations (50 nM - 6.25 nM) for 30s at 400 ml/min. Raw data were double referenced and globally fit with a 1:1 binding model showing accurate determination of an off-rate around  $10 \text{ s}^{-1}$ .


Sensorgrams of a 6-mer oligonucleotide (1.7 kDa) binding onto its complementary ssDNA (11 kDa biotinylated 34-mer) captured on streptavidin on a PCP-5 WAVEchip. The interaction was measured at 25°C. Zoom into the dissociation shows excellent data fitting and an accurate determination of an off-rate around  $10 \text{ s}^{-1}$ .

### SMALL MOLECULES CAN'T HIDE ANYMORE

Sensitivity is key and often limiting for accurate and reliable analysis of molecular interactions. The high-sensitivity of the Creoptix<sup>®</sup> WAVEsystem allows researchers to confidently analyze binding interactions at very low signal levels and high analyte-to-ligand molecular weight (MW) ratios.



Sensorgrams of acetazolamide (222 Da) binding to Carbonic Anhydrase II (29 kDa) immobilised at low density onto a PCH WAVEchip.



Sensorgrams of a small drug molecule (295 Da) binding to a target protein (110 kDa) immobilised at low density onto a PCH WAVEchip. Note the analyte-to-ligand MW ratio is  $> 300$ .

GCI IS FEATURED IN

**NATURE** Mechanistic basis for the activation of plant membrane receptor kinases by SERK-family coreceptors

**PLANT PHYSIOLOGY** The SERK3 elongated allele defines a role for BIR ectodomains in brassinosteroid signalling

**NATURE COMMUNICATIONS** Dynamics of human protein kinase Aurora A linked to drug selectivity

REFERENCES

1. Kozma et al., "Grating coupled optical waveguide interferometer for label-free biosensing", Sensors and Actuators B: Chemical, 155:444-450 (2011)
2. Höhmann et al., "The SERK3 elongated allele defines a role for BIR ectodomains in brassinosteroid signalling", Nature Plants, 4:345-351 (2018)
3. Brough et al., "Application of Off-Rate Screening in the Identification of Novel Pan-isoform Inhibitors of Pyruvate Dehydrogenase Kinase", J. Med. Chem., 40 (6), pp 2271-2285 (2017)
4. Brough et al., "4,5-Diaryloxazole Hsp90 Chaperone Inhibitors: Potential Therapeutic Agents for the Treatment of Cancer", J. Med. Chem., 51 (2), pp 190-218 (2008)

Scan this code to request a copy of this poster

