

Antibody Characterization in Physiologically Relevant Conditions by Grating-Coupled Interferometry (GCI)

Matyas Vegh¹, Jealemy Galindo¹, Fabio Andres¹, Fabio M. Spiga¹, Richard J. Isaacs²

1 Creoptix AG, 8820 Wädenswil, Switzerland, info@creoptix.com

2 Creoptix, Inc, 111 South Bedford Street, Ste.108, Boston, MA-01803, Burlington, USA



Abstract


Protein interaction analysis is crucial for the successful identification of drug targets. Amongst the wide range of biophysical techniques available today, surface-based biosensors such as SPR utilize optical systems to study biomolecular interactions without labels.

Grating-Coupled Interferometry (GCI)¹ is a novel label-free technique suitable for the analysis of low-molecular-mass drugs, multiprotein complexes and even larger particles like viruses or bacteria, with a remarkably wide affinity range (from mM to pM)².

Here, we show selected examples of GCI applications in the field of antibody characterization, including binding kinetics and affinity measurements in physiologically relevant media (serum, plasma, cell supernatant), that highlight how this technology is revolutionizing the study of molecular interactions and drug discovery.

Grating-Coupled Interferometry (GCI) technology

GCI is a surface-based, label-free biosensing technique. When target molecules (e.g. proteins γ) are attached to the sensor surface, binding analytes (\bullet) lead to an increase in mass and hence to a change in the refractive index within the evanescent field near the surface.



In GCI, refractive index changes on a sensor surface are measured as time-dependent phase-shift signals. The long-light to sample interaction length of the waveguide provides intrinsically high signal-to-noise levels³.



The Creoptix[®] WAVE system combines GCI with innovative no-clog microfluidics, allowing the study of interactions even between very small ligands and large receptors.


The system uses a robust sensor and microfluidics cartridge, the WAVEchip[®], where proteins, such as antibodies, VLPs, peptides or other molecules can be immobilized using various chemistries.



Kinetics in Serum and Plasma


Binding interaction analysis in physiologically relevant conditions is possible using GCI.

Monoclonal antibodies (mAbs) were immobilized via amine coupling on a 4PCH chip. The protein analyte was spiked in 90% serum or 90% plasma as indicated, and injected in six concentrations (4 nM - 31 pM). Raw data were double referenced and globally fit with a 1:1 binding model.



Kinetic Analysis on GPCRs

GCI can be used to kinetically characterize molecules binding to G-Protein-Coupled Receptors (GPCRs) directly from crude membrane preps or virus-like particles (VLPs).




CXCR4-VLPs (Integral Molecular) were captured on a 4PCP chip coated with Wheat Germ Agglutinin (WGA). Kinetics of a nanobody binding to its biotinylated GPCR target in a crude membrane prep (captured with WGA). Kinetics of an anti-CXCR4 mAb interacting with CXCR4-VLPs were determined by injecting five concentrations of the nanobody (0.96 nM to 3 μ M).

Antibody Profiling in Serum

The robust sensor and microfluidics of the technology allows the kinetic characterization of molecular interactions in almost pure serum and even plasma. Binding proteins developed for diagnostic applications can therefore be directly profiled in great depth within the relevant matrix.

Comprehensive interaction data directly acquired in early-stage antibody profiling enables valuable kinetic characterization of diagnostic leads in the physiologically relevant matrix. Here, two different antibodies were amine-coupled to a 4PCP chip and the antigen was injected in either buffer (PBS P⁺) or 90% human serum in a dilution series of eight concentrations (137 pM to 300 nM). Raw data were double referenced and globally fit with a 1:1 binding model.



Antibody 1 demonstrates its suitability for diagnostic application, revealing highly similar kinetics in 90% serum as in buffer. On the other hand, antibody 2 evidently shows no binding activity in the presence of serum and is therefore not suited for diagnostic application.

GCI is featured in

Mechanistic basis for the activation of plant membrane receptor kinases by SERK-family co-receptors

The SERK3 elongated allele defines a role for BIR ectodomains in brassinosteroid signalling

Dynamics of human protein kinase Aurora A linked to drug selectivity

References

- 1_P. Kozma et al., "Grating coupled optical waveguide interferometer for label-free biosensing", Sensors and Actuators B: Chemical, 155:446-450 (2011)
- 2_Hohmann et al., "The SERK3 elongated allele defines a role for BIR ectodomains in brassinosteroid signalling.", Nature Plants, 4:345-351 (2018)
- 3_www.sprpages.nl/how-10v

Anti-Drug Antibody Quantification from Blood Serum

Anti-drug antibody (ADA) titers can be quantified by GCI directly from blood serum for immunogenicity detection.

As shown here, in combination with robust microfluidics GCI allows for the detection of significantly lower antibody titers (lower limit of detection, LoD) by enabling direct quantification in almost pure blood serum.

