


---

# WAVEchips

October 7, 2015

---



# WAVEchips

All our sensor chips have the sensor surface sealed to the integrated microfluidics. This provides optimal sensor surface protection at most convenience.

Product Name	Description	Quantity	Order No
WAVEchip DXH	Versatile sensor chip with high capacity carboxymethyl-dextrane functionalization. Suited for immobilization with –NH <sub>2</sub> , –SH, –CHO, –OH or –COOH groups.	Pack of 1	DXH-01
		Pack of 3	DXH-03
		Pack of 10	DXH-10
WAVEchip DXP	Versatile sensor chip with low capacity carboxymethyl-dextrane functionalization. Suited for immobilization with –NH <sub>2</sub> , –SH, –CHO, –OH or –COOH groups.	Pack of 1	DXP-01
		Pack of 3	DXP-03
		Pack of 10	DXP-10
WAVEchip DXH-S	High capacity carboxymethyl-dextrane functionalization pre-immobilized with streptavidin for capturing biotinylated molecules	Pack of 1	DXH-S-01
		Pack of 3	DXH-S-03
		Pack of 10	DXH-S-10
WAVEchip PCH	Versatile sensor chip with high capacity polycarboxylate functionalization. Suited for immobilization with –NH <sub>2</sub> , –SH, –CHO, –OH or –COOH groups.	Pack of 1	PCH-01
		Pack of 3	PCH-03
		Pack of 10	PCH-10
WAVEchip PCL	Versatile sensor chip with low capacity polycarboxylate functionalization. Suited for immobilization with –NH <sub>2</sub> , –SH, –CHO, –OH or –COOH groups.	Pack of 1	PCL-01
		Pack of 3	PCL-03
		Pack of 10	PCL-10
WAVEchip PCP	Versatile sensor chip with a thin polycarboxylate functionalization. Suited for immobilization with –NH <sub>2</sub> , –SH, –CHO, –OH or –COOH groups. For molecules and larger particles.	Pack of 1	PCP-01
		Pack of 3	PCP-03
		Pack of 10	PCP-10

## WAVEchips, continued

Product Name	Description	Quantity	Order No
WAVEchip PCH-N	High capacity sensor chip with pre-immobilized NTA for capturing His tagged molecules	Pack of 1	PCH-N-01
		Pack of 3	PCH-N-03
		Pack of 10	PCH-N-10
WAVEchip PCH-S	High capacity sensor chip with pre-immobilized streptavidin for capturing biotinylated molecules	Pack of 1	PCH-S-01
		Pack of 3	PCH-S-03
		Pack of 10	PCH-S-10
WAVEchip PCP-S	Thin polycarboxylate sensor chip with pre-immobilized streptavidin for capturing biotinylated molecules or larger particles	Pack of 1	PCP-S-01
		Pack of 3	PCP-S-03
		Pack of 10	PCP-S-10